



# ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Steven R Fischer</i>		
Location <i>Ward County Jail</i>	Intoxilyzer® Serial Number <i>80-003066</i>	
Gas Lot Number <i>11611080A2</i>	Gas Cylinder Number <i>108</i>	Gas Expiration Date <i>06/01/2013</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.078*AC
  - B. *0.078*AC
  - C. *0.078*AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> #272, NDHP	Date <i>9-7-2011</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>02 Apr 2012</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003066  
Location = WARD                    8164.13.00 06/09  
09/07/2011                            14:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:37
02 Std. Gas	0.078	14:37
03 Room Air	0.000	14:38
04 Std. Gas	0.078	14:38
05 Room Air	0.000	14:39
06 Std. Gas	0.078	14:39
07 Room Air	0.000	14:40

Lot No = 11611080A2  
Cyl No = 108  
Exp Date = 06/01/2013  
County = 51                            Oper No. = 049904



Operator Signature  
STEVEN FISCHER

Remarks: