



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Jamie M Huschka</u>		
Location: <u>Williams Co LEC</u>		
Intoxilyzer® Serial Number: <u>80-004189</u>		
Gas Lot Number: <u>15410080A1</u>	Gas Cylinder Number: <u>6</u>	Expiration Date: <u>07/01/12</u>
<input checked="" type="checkbox"/> New Cylinder Installed <input type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: <u>0.079</u> %		
Test 2: <u>0.079</u> %		
Test 3: <u>0.078</u> %		
(ATTACH TEST RECORD)		

Jamie M Huschka
Chemical Test Operator Signature

04 Apr 11
Date

Deb Kashner
Reviewed By

13 Apr 11
Date

80-004189
P. 122 of 124
PM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = ~~TEXT~~ WILL 8164.13.00 06/09
04/04/2011 DK 13 Apr 11 17:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:54
02 Std. Gas	0.079	17:54
03 Room Air	0.000	17:55
04 Std. Gas	0.079	17:55
05 Room Air	0.000	17:55
06 Std. Gas	0.078	17:56
07 Room Air	0.000	17:56

Lot No = 15410080A1
Cyl No = 6
Exp Date = 07/01/2012
County = 53 Oper No. = 100406



Operator Signature
JAMIE M. HUSCHKA

Remarks:

Form 106-I8000

80-004189
P. 123 of 124
FM

13 April 2011

I emailed HP Jamie Huschka regarding 8000 Installation paperwork. He forgot to change the location code in the Intoxilyzer when he did an install 04 April 2011. I first called the office, but no one answered. After the email I called the office again to leave a message, but Darcy Aberle answered and he said he would go downstairs to fix this issue.

Deb Kashur

80-004189
p. 124 of 124
JFM