



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (1-2010)

Serial Number: 80-004938 Instrument Location: Watford City

Reason for Install/Repair:
 Install After Lab Repair/Inspection
 Other (Specify) _____

Check When Done:

- A. Surge Protector Installed/Properly Grounded.
- B. Telephone Line Connected to Intoxilyzer® 8000.
- C. Breath Tube Heated.
- D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- E. Review/Enter Gas Setup (Level 1, Function S).
- F. Print Test (Level 1, Function P). (Attach Test Record.)
- G. ACA Test (Level 1, Function C). (Attach Test Record.)
- H. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- I. Repair and/or Maintenance Performed: _____

Deborah N...
 Field Inspector Signature

4/6/11
 Date

Deb Kashur
 Reviewed By

14 Apr 11
 Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004938
Location = WATF 8164.13.00 06/09
04/06/2011 22:58

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:59
02 Std. Gas	0.080	22:59
03 Room Air	0.000	23:00
04 Std. Gas	0.080	23:00
05 Room Air	0.000	23:01
06 Std. Gas	0.081	23:01
07 Room Air	0.000	23:02

Lot No = 659357
Cyl No = 44
Exp Date = 09/03/2011
County = 27 Oper No. = 030924



Operator Signature
DUSTIN B. NEWMAN

Remarks:

Form 106-I8000

80-004938
P. 84 of 85
JEM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004938
Location = WATF 8164.13.00 06/09
04/06/2011 23:07

Test	AC	Time
01 Diagnostic	OK	23:11
02 Room Air	RFI*	23:12
03 Room Air	RFI*	23:12

*Invalid Test
Inhibited - RFI

Sub Name = TEST, SUBJECT NA
Sub DOB = 03/25/1980
Sub Sex = Male Weight = 215
Test = OTH Cit = NA
Dr. Lic. = ND/NA
Lot No = 659357
Cyl No = 44
Expiration Date = 09/03/2011
County = 27 Oper No. = 030924

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DUSTIN B. NEWMAN

Remarks: RFI Test

Form 106-I8000

80-004938
p. 85 of 85
JMN