

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004204

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN MP3066 Lot No. 15160 Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3068 Lot No. 14102 Rep AC 0.400 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. ICS 4 Solution 0.10AC+0.05% ACETONE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 13314080A3 Cyl No. 028 Exp. Date 6/5/2016

Test 1 <u>0.078</u>	Test 4 <u>0.078</u>	Test 7 <u>0.078</u>
Test 2 <u>0.078</u>	Test 5 <u>0.079</u>	Test 8 <u>0.078</u>
Test 3 <u>0.079</u>	Test 6 <u>0.078</u>	Test 9 <u>0.079</u>

Average 0.078

J. Tank pressure, Level 3 Func D Sub G P: 244 psi Reg ~250 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION, REPLACED

SIMULATED RETURN O-RING

4.14.2016  
Date

14 Apr 2016  
Date

Charles E. Eden  
Field Inspector's Signature

Deb Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                 8164.13.00 06/09  
04/14/2016                       10:22

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:                Enabled  
Start Test Sequence:            DABACABA  
Display Prelim Rslt?            Yes  
Display Third Digit?            Yes  
Inhib Printer(Y/N)?            No  
Display Volume?                 No  
Disable On Memfull?            Yes  
# of Print Copies?              1  
Select Std (D/W/I)?            Dry  
Standard Value?                 0.080  
Standard Lot #?                 13314080A3  
Standard Cyl #?                 28  
Standard Expiration?            06/05/2016  
Oper No?                         666666

Flow Cal. Date:                 07/13/2010  
Slope                            750  
Intercept                       -815939

IR Calibration Date:            05/12/2015

	3um	9um
0th Coef(*100):	-23668	-20605
1st Coef(*100):	259586	131325
2nd Coef(*100):	3023	1492
H2O adj(mg/l*10k):	445	418

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                 8164.13.00 06/09  
04/14/2016                        10:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:25
02 Std. Sol.	0.018	10:26
03 Room Air	0.000	10:26
04 Std. Sol.	0.019	10:27
05 Room Air	0.000	10:28
06 Std. Sol.	0.019	10:28
07 Room Air	0.000	10:29

08 Sim Temp = 34.0°C

Simul Ser No = MP3066

Std Sol No = 15160

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC 0.0204C

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004204  
Location = TOXL      8164.13.00 06/09  
04/14/2016      10:31

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:31
02 Std. Sol.	0.396	10:32
03 Room Air	0.000	10:33
04 Std. Sol.	0.395	10:33
05 Room Air	0.000	10:34
06 Std. Sol.	0.396	10:35
07 Room Air	0.000	10:35

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 14102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGH AC*  
*0.400 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                 8164.13.00 06/09  
04/14/2016                        10:41

Test	AC	Time
01 Room Air	0.000	10:43
02 *Subject Test	INT*	10:44
03 Room Air	0.000	10:45

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                 Weight = 150  
Test = DUI                         Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                8164.13.00 06/09  
04/14/2016                      10:45

Test	AC	Time
01 Diagnostic	OK	10:46
02 Room Air	0.000	10:46
03 *Subject Test	RFI*	10:47
04 Room Air	0.000	10:47

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                Weight = 150  
Test = DUI                      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08                      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                    8164.13.00 06/09  
04/14/2016                            10:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:48
02 Std. Gas	0.078	10:48
03 Room Air	0.000	10:49
04 Std. Gas	0.078	10:49
05 Room Air	0.000	10:50
06 Std. Gas	0.079	10:50
07 Room Air	0.000	10:51

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                    8164.13.00 06/09  
04/14/2016                            10:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:51
02 Std. Gas	0.078	10:52
03 Room Air	0.000	10:52
04 Std. Gas	0.079	10:53
05 Room Air	0.000	10:53
06 Std. Gas	0.078	10:54
07 Room Air	0.000	10:54

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004204  
Location = TOXL                8164.13.00 06/09  
04/14/2016                      10:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:55
02 Std. Gas	0.078	10:55
03 Room Air	0.000	10:56
04 Std. Gas	0.078	10:56
05 Room Air	0.000	10:57
06 Std. Gas	0.079	10:57
07 Room Air	0.000	10:58

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08                      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC