



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>David Feland</u>		
Location: <u>Kidder SO.</u>		
Intoxilyzer® Serial Number: <u>80-004186</u>		
Gas Lot Number: <u>659357</u>	Gas Cylinder Number: <u>37</u>	Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: <u>0.081</u> %		
Test 2: <u>0.080</u> %		
Test 3: <u>0.080</u> %		
(ATTACH TEST RECORD)		

David Feland
Chemical Test Operator Signature

29 Apr 11
Date

Deb Kashner
Reviewed By

06 May 2011
Date

80-004186
P. 104 of 105
km

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = STEE 8164.13.00 06/09
04/29/2011 15:01

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:01
02 Std. Gas	0.081	15:02
03 Room Air	0.000	15:02
04 Std. Gas	0.080	15:03
05 Room Air	0.000	15:03
06 Std. Gas	0.080	15:04
07 Room Air	0.000	15:04

Lot No = 659357
Cyl No = 37
Exp Date = 09/03/2011
County = 22 Oper No. = 119404


Operator Signature
DAVID FELAND

Remarks:

Form 106-I8000

80-004186
P.105 of 105
JFM