



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>CHARLES E. EDER</i>		
Location <i>ASHL</i>	Intoxilyzer® Serial Number <i>80-005360</i>	
Gas Lot Number <i>3391308044</i>	Gas Cylinder Number <i>041</i>	Gas Expiration Date <i>01/01/2016</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.081* AC
 - B. *0.083* AC
 - C. *0.082* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Charles E. Eder</i>	Date <i>02 MAY 14</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>05 May 14</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/02/2014 10:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:49
02 Std. Gas	0.081	10:49
03 Room Air	0.000	10:50
04 Std. Gas	0.083	10:50
05 Room Air	0.000	10:51
06 Std. Gas	0.082	10:51
07 Room Air	0.000	10:51

Lot No = 33913080A4

Cyl No = 41

Exp Date = 01/01/2016

County = 26

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CYLINDER INSTALL
0.080 AC

Form 106-I8000