

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-003066

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5132 Lot No. 201507F Rep AC 0.020 ✓
- F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN DR3841 Lot No. 201603A Rep AC 0.300 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5144 Lot No. ICS4 Solution 0.10 AC + 0.05% acetone ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 14214080A6 Cyl No. 51 Exp. Date 07/05/2016

Test 1 0.079 Test 4 0.080 Test 7 0.079

Test 2 0.080 Test 5 0.079 Test 8 0.079

Test 3 0.079 Test 6 0.079 Test 9 0.080

Average 0.079

J. Tank pressure, Level 3 Func D Sub G P: 509 psi Reg 500 psi

K. Remarks/Maintenance record: Annual inspection

29 Apr 2016  
Date

02 MAY 2016  
Date

Deb Shanover  
Field Inspector's Signature

Charles Elder  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003066  
Location = TOXL                    8164.13.00 06/09  
04/29/2016                            11:38

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        No  
# of Print Copies?           1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            14214080A6  
Standard Cyl #?            51  
Standard Expiration?        07/05/2016  
Oper No?                    777777

Flow Cal. Date:            07/29/2010  
Slope                        806  
Intercept                  -895033

IR Calibration Date:        07/29/2010  
                                  3um                    9um

-----  
0th Coef(\*100):            -14075                -5808  
1st Coef(\*100):            272113                138204  
2nd Coef(\*100):            1637                   986  
H2O adj(mg/l\*10k):        289                    148

\*\*\*\*\* Printer Test End \*\*\*\*\*

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Print test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      15:46

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:47
02 Std. Sol.	0.021	15:48
03 Room Air	0.000	15:48
04 Std. Sol.	0.021	15:49
05 Room Air	0.000	15:49
06 Std. Sol.	0.021	15:50
07 Room Air	0.000	15:51

08 Sim Temp = 34.0°C

Simul Ser No = DR5132

Std Sol No = 201507F

County = 08

Oper No. = 777777



Operator Signature

DEB SHANAVER

Remarks:

Low AC      0.021 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      15:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:42
02 Std. Sol.	0.298	15:42
03 Room Air	0.000	15:43
04 Std. Sol.	0.300	15:44
05 Room Air	0.000	15:44
06 Std. Sol.	0.300	15:45
07 Room Air	0.000	15:45

08 Sim Temp = 34.0°C

Simul Ser No = DR3841  
Std Sol No = 201603A  
County = 08      Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*High Ac      0.300AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      15:54

Test	AC	Time
01 Diagnostic	OK	15:55
02 Room Air	RFI*	15:55
03 Room Air	0.000	15:55

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 14214080A6  
Cyl No = 51  
Expiration Date = 07/05/2016  
County = 08      Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
DEB SHANAVER

Remarks:

*RFI test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      11:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:35
02 Std. Gas	0.079	11:36
03 Room Air	0.000	11:36
04 Std. Gas	0.080	11:36
05 Room Air	0.000	11:37
06 Std. Gas	0.079	11:37
07 Room Air	0.000	11:38

Lot No = 14214080A6  
Cyl No = 51  
Exp Date = 07/05/2016  
County = 08      Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Calibration check 0.080AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      11:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:39
02 Std. Gas	0.080	11:39
03 Room Air	0.000	11:40
04 Std. Gas	0.079	11:40
05 Room Air	0.000	11:41
06 Std. Gas	0.079	11:41
07 Room Air	0.000	11:42

Lot No = 14214080A6  
Cyl No = 51  
Exp Date = 07/05/2016  
County = 08      Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Calibration check 0.080AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = TOXL      8164.13.00 06/09  
04/29/2016      15:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:30
02 Std. Gas	0.079	15:31
03 Room Air	0.000	15:31
04 Std. Gas	0.079	15:32
05 Room Air	0.000	15:32
06 Std. Gas	0.080	15:32
07 Room Air	0.000	15:33

Lot No = 14214080A6  
Cyl No = 51  
Exp Date = 07/05/2016  
County = 08      Oper No. = 777777

*Deb Shanaver*

Operator Signature  
DEB SHANAVER

Remarks:

*Calibration check 0.080AC*

Form 106-I8000