

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003058

Location: Ward Co

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u> </u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5115</u> Lot No. <u> </u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u> </u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |

K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. <u>459358</u>	Cyl No. <u>32</u>	Exp. Date <u>09/03/11</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.082</u>
Test 2 <u>0.082</u>	Test 5 <u>0.081</u>	Test 8 <u>0.082</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.082</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 645 psi Reg 645 psi

M. Remarks/Maintenance record: Annual inspection

06 May 2010
Date

Deb Kasher
Field Inspector's Signature

15 May 2010
Date

Mary Plawn
Reviewed by

80-003058
P. 84492
JK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:33

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 32
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 04/07/2008
Slope 698
Intercept -678970

IR Calibration Date: 06/05/2008
 3um 9um

0th Coef(*100):	-18867	-13178
1st Coef(*100):	274682	135660
2nd Coef(*100):	1139	678
H2O adj(mg/l*10k):	434	216

***** Printer Test End *****


Operator Signature
DEB KASHUR

Remarks:

Print Test

Form 106-I8000

80-003058
P. 85 of 92
YK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:50

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:51
02 Std. Sol.	0.020	09:52
03 Room Air	0.000	09:52
04 Std. Sol.	0.020	09:53
05 Room Air	0.000	09:53
06 Std. Sol.	0.020	09:54
07 Room Air	0.000	09:55

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 51

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021AC

Form 106-I8000

80-003058
P. 86 of 92
JKM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 10:00

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:01
02 Std. Sol.	0.288	10:01
03 Room Air	0.000	10:02
04 Std. Sol.	0.290	10:03
05 Room Air	0.000	10:03
06 Std. Sol.	0.291	10:04
07 Room Air	0.000	10:04

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 51

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

DK
Calibration Linearity test 0.298 AC

Form 106-I8000

80-003058
P. 87-192
JFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:46

Test	AC	Time
01 Room Air	0.000	09:47
02 *Subject Test	INT*	09:47
03 Room Air	0.000	09:48

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 32
Expiration Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Interferent check

Form 106-I8000

*0.00 AC
+acetone*

80-003058
P. 88 of 92
FA

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:48

Test	AC	Time
01 Diagnostic	OK	09:49
02 Room Air	RFI*	09:49
03 Room Air	0.000	09:50

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 32
Expiration Date = 09/03/2011
County = 51 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

80-003058
P. 89/92
JCM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:35
02 Std. Gas	0.081	09:36
03 Room Air	0.000	09:36
04 Std. Gas	0.082	09:37
05 Room Air	0.000	09:37
06 Std. Gas	0.081	09:37
07 Room Air	0.000	09:38

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003058
P. 90492
KH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:39
02 Std. Gas	0.081	09:40
03 Room Air	0.000	09:40
04 Std. Gas	0.081	09:40
05 Room Air	0.000	09:41
06 Std. Gas	0.081	09:41
07 Room Air	0.000	09:42

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003058
P. 9/1/92
KSM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003058
Location = WARD 8164.13.00 06/09
05/06/2010 09:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:42
02 Std. Gas	0.082	09:43
03 Room Air	0.000	09:43
04 Std. Gas	0.082	09:44
05 Room Air	0.000	09:44
06 Std. Gas	0.082	09:44
07 Room Air	0.000	09:45

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003058
p. 92 of 92
YFM