

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004198

Location: Renoville Co - Mohale

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>---</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5115</u> Lot No. <u>---</u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u>---</u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659358</u>	Cyl No. <u>27</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.079</u>	Test 5 <u>0.079</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.078</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 825 psi Reg 820 psi

M. Remarks/Maintenance record: Annual inspection

06 May 2010
Date

Deb Kashur
Field Inspector's Signature

15 May 2010
Date

Mary Pearson
Reviewed by

80-004198
P. 84/92 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:31

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 27
Standard Expiration? 09/03/2011
Oper No? 777777

Flow Cal. Date: 06/04/2009
Slope 671
Intercept -536290

IR Calibration Date: 12/09/2009
3um 9um

	3um	9um
0th Coef(*100):	-31070	-34849
1st Coef(*100):	268664	135239
2nd Coef(*100):	-643	440
H2O adj(mg/l*10k):	570	615

***** Printer Test End *****

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-I8000

80-004198
P. 85/92 VK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:48
02 Std. Sol.	0.019	12:49
03 Room Air	0.000	12:49
04 Std. Sol.	0.018	12:50
05 Room Air	0.000	12:50
06 Std. Sol.	0.018	12:51
07 Room Air	0.000	12:52

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 38

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021 AC

Form 106-18000

80-004198
P. 8/2 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 13:00

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:01
02 Std. Sol.	0.285	13:02
03 Room Air	0.000	13:02
04 Std. Sol.	0.286	13:03
05 Room Air	0.000	13:03
06 Std. Sol.	0.285	13:04
07 Room Air	0.000	13:05

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 38

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.298AC

Form 106-I8000

80-004198
P. 87/92 *JK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:42

Test	AC	Time
01 Room Air	0.000	12:43
02 *Subject Test	INT*	12:44
03 Room Air	0.000	12:44

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 27
Expiration Date = 09/03/2011
County = 38 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Interferent check

*0.100 AC
+ acetone*

Form 106-18000

80-004198
P. 88/92 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:44

Test	AC	Time
01 Diagnostic	OK	12:45
02 Room Air	RFI*	12:45
03 Room Air	0.000	12:46

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 27
Expiration Date = 09/03/2011
County = 38 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

80-004198
P. 89/92 *AK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:33
02 Std. Gas	0.079	12:33
03 Room Air	0.000	12:34
04 Std. Gas	0.079	12:34
05 Room Air	0.000	12:34
06 Std. Gas	0.079	12:35
07 Room Air	0.000	12:35

Lot No = 659358

Cyl No = 27

Exp Date = 09/03/2011

County = 38

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

*80-004198
P. 90/92 YH*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:36
02 Std. Gas	0.079	12:36
03 Room Air	0.000	12:37
04 Std. Gas	0.079	12:37
05 Room Air	0.000	12:38
06 Std. Gas	0.078	12:38
07 Room Air	0.000	12:39

Lot No = 659358

Cyl No = 27

Exp Date = 09/03/2011

County = 38

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004198
P. 9/1/92 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = MOHL 8164.13.00 06/09
05/06/2010 12:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:39
02 Std. Gas	0.079	12:40
03 Room Air	0.000	12:40
04 Std. Gas	0.079	12:41
05 Room Air	0.000	12:41
06 Std. Gas	0.079	12:41
07 Room Air	0.000	12:42

Lot No = 659358

Cyl No = 27

Exp Date = 09/03/2011

County = 38

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

*80-004198
P. 92/92 KHL*