



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Richard A. Sherlock (4713)</u>		
Location: <u>GRAF</u>		
Intoxilyzer® Serial Number: <u>80-003072</u>		
Gas Lot Number: <u>659358</u>	Gas Cylinder Number: <u>16</u>	Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>081</u> %		
Test 2: 0. <u>082</u> %		
Test 3: 0. <u>081</u> %		
(ATTACH TEST RECORD)		

Richard A. Sherlock (4713)
 Chemical Test Operator Signature

05/08/10
 Date

Deb Kashur
 Reviewed By

11 May 2010
 Date

80-003072
 7/29/10


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003072
Location = GRAF 8164.13.00 06/09
05/08/2010 15:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:36
02 Std. Gas	0.081	15:36
03 Room Air	0.000	15:37
04 Std. Gas	0.082	15:37
05 Room Air	0.000	15:38
06 Std. Gas	0.081	15:38
07 Room Air	0.000	15:39

Lot No = 659358
Cyl No = 16
Exp Date = 09/03/2011
County = 50 Oper No. = 100428



Operator Signature
RICHARD A. SHERLOCK

Remarks: *ACA*

Form 106-I8000

80-003072
79/79 *KAM*