



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

| | | |
|--|---------------------------|---|
| Chemical Test Operator Name (Print) Daniel C. Garvens | | |
| Location Pembina County Jail (CAVA) | | Intoxilyzer® Serial Number 80-005943 |
| Gas Lot Number 13314080A3 | Gas Cylinder Number 33 | Gas Expiration Date 06/05/2016 |

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.029 AC
 - B. 0.029 AC
 - C. 0.029 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| | |
|--|-------------------|
| Chemical Test Operator Signature <i>[Signature]</i> | Date 5-8-15 |
| Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i> | Date 14 May 15 |

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005943
Location = CAVA 8164.13.00 06/09
05/08/2015 16:36

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:36 |
| 02 Std. Gas | 0.079 | 16:37 |
| 03 Room Air | 0.000 | 16:37 |
| 04 Std. Gas | 0.079 | 16:38 |
| 05 Room Air | 0.000 | 16:38 |
| 06 Std. Gas | 0.079 | 16:39 |
| 07 Room Air | 0.000 | 16:39 |

Lot No = 13314080A3
Cyl No = 33
Exp Date = 06/05/2016
County = 34 Oper No. = 060849



Operator Signature
DANIEL C. GARVENS

Remarks: *New Intox install*