



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number 80-004196	Instrument Location MINOT AFB
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Malthe [Signature]</i>	Date 10 MAY 2013
Reviewed By (Crime Laboratory Use Only) <i>Deb Keshur</i>	Date 14 May 2013

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004196
Location = MNAF 8164.13.00 06/09
05/10/2013 17:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:22
02 Std. Gas	0.078	17:22
03 Room Air	0.000	17:23
04 Std. Gas	0.078	17:23
05 Room Air	0.000	17:24
06 Std. Gas	0.079	17:24
07 Room Air	0.000	17:24

Lot No = 16512080A2
Cyl No = 38
Exp Date = 08/01/2014
County = 64 Oper No. = 129952



Operator Signature
MATTHEW STASO

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004196
Location = MNAF 8164.13.00 06/09
05/10/2013 17:31

Test	AC	Time
01 Diagnostic	OK	17:33
02 Room Air	RFI*	17:33
03 Room Air	RFI*	17:33

*Invalid Test
Inhibited - RFI

Sub Name = BROOKS III, ROBERT LEWIS
Sub DOB = 09/09/1990
Sub Sex = Male Weight = 185
Test = OTH Cit = NA
Dr. Lic. = ND/BRO903384
Lot No = 16512080A2
Cyl No = 38
Expiration Date = 08/01/2014
County = 64 Oper No. = 129952

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MATTHEW STASO

Remarks:

Form 106-I8000