

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 56
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 06/23/2009
Slope 681
Intercept -622430

IR Calibration Date: 12/10/2009
 3um 9um

0th Coef(*100): -30896 -34419
1st Coef(*100): 274294 133969
2nd Coef(*100): -715 552
H2O adj(mg/l*10k): 491 645

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:30
02 Std. Sol.	0.016	10:31
03 Room Air	0.000	10:32
04 Std. Sol.	0.016	10:32
05 Room Air	0.000	10:33
06 Std. Sol.	0.016	10:34
07 Room Air	0.000	10:34

08 Sim Temp = 34.0°C

Simul Ser No = 000000000
Std Sol No = 00000
County = 62 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low 0.021 AC

Form 106-18000

*80-004194
P. 80/86 FR*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:36

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:36
02 Std. Sol.	0.287	10:37
03 Room Air	0.000	10:38
04 Std. Sol.	0.288	10:39
05 Room Air	0.000	10:39
06 Std. Sol.	0.289	10:40
07 Room Air	0.000	10:41


08 Sim Temp = 34.0°C

Simul Ser No = 000000000

Std Sol No = 00000

County = 62

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-I8000

80-004194
P. 8/86 JCM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:17

Test	AC	Time
01 Room Air	0.000	10:18
02 *Subject Test	INT*	10:18
03 Room Air	0.000	10:19

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 000000000
Std Sol No = 00000
County = 62 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Interferent 0.1ACT 0.05% Acetone

Form 106-I8000

*80-004194
P. 8/86 KM*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:13

Test	AC	Time
01 Diagnostic	OK	10:14
02 Room Air	RFI*	10:15
03 Room Air	0.000	10:15

*Invalid Test
Inhibited - RFI

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004194
P. 83/86 JCH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:03
02 Std. Gas	0.078	10:03
03 Room Air	0.000	10:04
04 Std. Gas	0.078	10:04
05 Room Air	0.000	10:04
06 Std. Gas	0.078	10:05
07 Room Air	0.000	10:05

Lot No = 659358
Cyl No = 56
Exp Date = 09/03/2011
County = 62 Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004194
P. 84/86 KA


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:07
02 Std. Gas	0.077	10:07
03 Room Air	0.000	10:07
04 Std. Gas	0.077	10:08
05 Room Air	0.000	10:08
06 Std. Gas	0.077	10:09
07 Room Air	0.000	10:09

Lot No = 659358
Cyl No = 56
Exp Date = 09/03/2011
County = 62 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

*80-004194
P. 85/86
KM*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = FTOT 8164.13.00 06/09
05/13/2010 10:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:43
02 Std. Gas	0.079	10:44
03 Room Air	0.000	10:44
04 Std. Gas	0.079	10:44
05 Room Air	0.000	10:45
06 Std. Gas	0.079	10:45
07 Room Air	0.000	10:46

Lot No = 659358
Cyl No = 56
Exp Date = 09/03/2011
County = 62 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

*80-004194
P. 8/26/10*