

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004952

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u> </u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5114</u> Lot No. <u> </u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u> </u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>709352</u>	Cyl No. <u>22</u>	Exp. Date <u>12/23/2011</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.079</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.079</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 633 psi Reg 640 psi

M. Remarks/Maintenance record: Annual inspection

14 May 2010
Date

Deb Kashur
Field Inspector's Signature

15 May 2010
Date

Mary Plaron
Reviewed by

80-004952

P. 8/16 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/14/2010 08:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:42
02 Std. Sol.	0.017	08:43
03 Room Air	0.000	08:43
04 Std. Sol.	0.018	08:44
05 Room Air	0.000	08:44
06 Std. Sol.	0.018	08:45
07 Room Air	0.000	08:46

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021 AC

Form 106-I8000

*80-004952
P. 10/16 YJM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/14/2010 08:48

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:49
02 Std. Sol.	0.289	08:50
03 Room Air	0.000	08:50
04 Std. Sol.	0.291	08:51
05 Room Air	0.000	08:52
06 Std. Sol.	0.292	08:52
07 Room Air	0.000	08:53

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

linearity test 0.298 AC

Form 106-18000

80-004952
P. 1/16 192

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/14/2010 08:54

Test	AC	Time
01 Room Air	0.000	08:55
02 *Subject Test	INT*	08:55
03 Room Air	0.000	08:56

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Simul Ser No = 00000
Std Sol No = 000
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.100 AC + acetone

Form 106-18000

80-004952
P. 12/16 PM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/13/2010 16:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:04
02 Std. Gas	0.079	16:05
03 Room Air	0.000	16:05
04 Std. Gas	0.079	16:06
05 Room Air	0.000	16:06
06 Std. Gas	0.079	16:06
07 Room Air	0.000	16:07

Lot No = 709352

Cyl No = 22

Exp Date = 12/23/2011

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004952
P. 14/16 VM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/13/2010 16:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:20
02 Std. Gas	0.079	16:21
03 Room Air	0.000	16:21
04 Std. Gas	0.080	16:21
05 Room Air	0.000	16:22
06 Std. Gas	0.080	16:22
07 Room Air	0.000	16:23

Lot No = 709352
Cyl No = 22
Exp Date = 12/23/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004952
P. 15/16 AK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = TOXL 8164.13.00 06/09
05/13/2010 16:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:26
02 Std. Gas	0.080	16:27
03 Room Air	0.000	16:27
04 Std. Gas	0.080	16:27
05 Room Air	0.000	16:28
06 Std. Gas	0.081	16:28
07 Room Air	0.000	16:29

Lot No = 709352

Cyl No = 22

Exp Date = 12/23/2011

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-18000

80-004952
P. 14/16 JPM