

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004946

Location: Emmons Co SO - LINT

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>12016</u> Rep AC <u>0.020</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3803</u> Lot No. <u>12050</u> Rep AC <u>0.300</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR384</u> Lot No. <u> </u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>23411080A1</u>	Cyl No. <u>9</u>	Exp. Date <u>10/01/2013</u>
Test 1 <u>0.081</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 1005 psi Reg 1000 psi

M. Remarks/Maintenance record: Annual inspection

18 May 2012
Date

Debi Kashur
Field Inspector's Signature

18 May 12
Date

Charles E. Egan
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:35

***** Printer Test *****

abcdefghijklmnopqrstu vwxyz1234567890-=_|
ABCDEFGHIJKLMN O PQRSTU VWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstu vwxyz1234567890-=_|
ABCDEFGHIJKLMN O PQRSTU VWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 9
Standard Expiration? 10/01/2013
Oper No? 777777

Flow Cal. Date: 01/18/2010
Slope 717
Intercept -728279

IR Calibration Date: 05/19/2011
 3um 9um

	3um	9um
0th Coef(*100):	-2231	-4869
1st Coef(*100):	282398	143476
2nd Coef(*100):	146	780
H2O adj(mg/l*10k):	303	375

***** Printer Test End *****

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:48
02 Std. Sol.	0.022	09:48
03 Room Air	0.000	09:49
04 Std. Sol.	0.022	09:50
05 Room Air	0.000	09:50
06 Std. Sol.	0.023	09:51
07 Room Air	0.000	09:51

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 15

Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

low AC 0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:52

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:53
02 Std. Sol.	0.307	09:54
03 Room Air	0.000	09:55
04 Std. Sol.	0.309	09:56
05 Room Air	0.000	09:56
06 Std. Sol.	0.309	09:57
07 Room Air	0.000	09:57

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 15

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:59

Test	AC	Time
01 Room Air	0.000	09:59
02 *Subject Test	INT*	10:00
03 Room Air	0.000	10:00

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 9
Expiration Date = 10/01/2013
County = 15 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.x0Ac + acetone

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 10:01

Test	AC	Time
01 Diagnostic	OK	10:02
02 Room Air	RFI*	10:02
03 Room Air	0.000	10:03

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 23411080A1
Cyl No = 9
Expiration Date = 10/01/2013
County = 15 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:36
02 Std. Gas	0.081	09:36
03 Room Air	0.000	09:37
04 Std. Gas	0.080	09:37
05 Room Air	0.000	09:38
06 Std. Gas	0.080	09:38
07 Room Air	0.000	09:39

Lot No = 23411080A1

Cyl No = 9

Exp Date = 10/01/2013

County = 15

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:39
02 Std. Gas	0.080	09:40
03 Room Air	0.000	09:40
04 Std. Gas	0.081	09:41
05 Room Air	0.000	09:41
06 Std. Gas	0.080	09:41
07 Room Air	0.000	09:42

Lot No = 23411080A1
Cyl No = 9
Exp Date = 10/01/2013
County = 15 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:
Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004946
Location = LINT 8164.13.00 06/09
05/18/2012 09:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:43
02 Std. Gas	0.081	09:43
03 Room Air	0.000	09:44
04 Std. Gas	0.080	09:44
05 Room Air	0.000	09:45
06 Std. Gas	0.080	09:45
07 Room Air	0.000	09:46

Lot No = 23411080A1
Cyl No = 9
Exp Date = 10/01/2013
County = 15 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC