



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 50281 (1-2010)

Serial Number: <u>80-604204</u>	Instrument Location: <u>Adams Co So HETT</u>
Reason for Install/Repair:	
<input checked="" type="checkbox"/> Install After Lab Repair/Inspection <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- A. Surge Protector Installed/Properly Grounded.
- B. Telephone Line Connected to Intoxilyzer® 8000.
- C. Breath Tube Heated.
- D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- E. Review/Enter Gas Setup (Level 1, Function S).
- F. Print Test (Level 1, Function P). (Attach Test Record.)
- G. ACA Test (Level 1, Function C). (Attach Test Record.)
- H. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- I. Repair and/or Maintenance Performed: _____

J M Mollat
 Field Inspector Signature

5-20-2011 DK
5-12-2011
 Date

Deb Keshu
 Reviewed By

23 May 2011
 Date



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (1-2010)

Serial Number: 80-604204 Instrument Location: Adams Co. SO HETT

Reason for Install/Repair:

Install After Lab Repair/Inspection

Other (Specify) _____

Check When Done:

- A. Surge Protector Installed/Properly Grounded.
- B. Telephone Line Connected to Intoxilyzer® 8000.
- C. Breath Tube Heated.
- D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- E. Review/Enter Gas Setup (Level 1, Function S).
- F. Print Test (Level 1, Function P). (Attach Test Record.)
- G. ACA Test (Level 1, Function C). (Attach Test Record.)
- H. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- I. Repair and/or Maintenance Performed: _____

 Field Inspector Signature Date

 Reviewed By Date

80-004204
 P. 106-111
 RM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004204
Location = HETT 8164.13.00 06/09
05/12/2011 14:57

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 33
Standard Expiration? 07/01/2012
Oper No? 049130

Flow Cal. Date: 07/13/2010
Slope 750
Intercept -815939

IR Calibration Date: 07/13/2010

	3um	9um
0th Coef(*100):	-24262	-22939
1st Coef(*100):	266889	136772
2nd Coef(*100):	1535	960
H2O adj(mg/l*10k):	459	425

***** Printer Test End *****


Operator Signature
EUGENE M. MOLBERT

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004204
Location = HETT 8164.13.00 06/09
05/12/2011 14:57

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 33
Standard Expiration? 07/01/2012
Oper No? 049130

Flow Cal. Date: 07/13/2010
Slope 750
Intercept -815939

IR Calibration Date: 07/13/2010

	3um	9um
0th Coef(*100):	-24262	-22939
1st Coef(*100):	266889	136772
2nd Coef(*100):	1535	960
H2O adj(mg/l*10k):	459	425

***** Printer Test End *****

Operator Signature
EUGENE M. MOLBERT

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004204
Location = HETT 8164.13.00 06/09
05/20/2011 19:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:48
02 Std. Gas	0.079	19:49
03 Room Air	0.000	19:49
04 Std. Gas	0.079	19:49
05 Room Air	0.000	19:50
06 Std. Gas	0.079	19:50
07 Room Air	0.000	19:51

Lot No = 15410080A1
Cyl No = 33
Exp Date = 07/01/2012
County = 01 Oper No. = 049130


Operator Signature
EUGENE M. MOLBERT

Remarks:

Form 106-I8000

80-004204
p.109 of 111
KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004204
Location = HETT 8164.13.00 06/09
05/12/2011 15:04

Test	AC	Time
01 Room Air	RFI*	15:05
02 Room Air	RFI*	15:05

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 07/07/1952
Sub Sex = Female Weight = 180
Test = DUI Cit = 57657657
Dr. Lic. = ND/7476476
Lot No = 15410080A1
Cyl No = 33
Expiration Date = 07/01/2012
County = 01 Oper No. = 049130



Operator Signature
EUGENE M. MOLBERT

Remarks:

Form 106-18000

80-004204
P.1104/111
JEM

