



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Kyle Johnson 130781</i>		
Location <i>800 4th Ave. E. West Fargo, ND</i>	Intoxilyzer® Serial Number <i>80-004944</i>	
Gas Lot Number <i>01316080A2</i>	Gas Cylinder Number <i>35</i>	Gas Expiration Date <i>02/05/2018</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. *0.081* AC
- B. *0.081* AC
- C. *0.081* AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature] #664</i>	Date <i>05/22/2016</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date <i>27 May 16</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = WFPD 8164.13.00 06/09
05/22/2016 13:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:22
02 Std. Gas	0.081	13:22
03 Room Air	0.000	13:23
04 Std. Gas	0.081	13:23
05 Room Air	0.000	13:24
06 Std. Gas	0.081	13:24
07 Room Air	0.000	13:25

Lot No = 01316080A2
Cyl No = 35
Exp Date = 02/05/2018
County = 09 Oper No. = 130781

 #664
Operator Signature
KYLE JOHNSON

Remarks: