

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004191

Location: TOXL

- | | Check When Done |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| A. Is the breath tube heated? | <u>✓</u> |
| B. Is the diagnostic check complete? | <u>✓</u> |
| C. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| D. Print test. (Attach test record.) | <u>✓</u> |
| E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3374</u> Lot No. <u>201501B</u> Rep AC <u>0.020</u> | <u>✓</u> |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3841</u> Lot No. <u>201603A</u> Rep AC <u>0.300</u> | <u>✓</u> |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u>ICS4</u> Solution <u>0.10 AC + 0.05% acetone</u> | <u>✓</u> |
| H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) | <u>✓</u> |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) <u>0.080</u> AC | |

Lot No. 14214080A16 Cyl No. 51 Exp. Date 07/05/2016

Test 1 0.079 Test 4 0.080 Test 7 0.079

Test 2 0.079 Test 5 0.079 Test 8 0.080

Test 3 0.080 Test 6 0.080 Test 9 0.079

Average 0.079

J. Tank pressure, Level 3 Func D Sub G P: 401 psi Reg ~420 psi

K. Remarks/Maintenance record: Annual inspection, replaced breath hose/simulator return o-rings, cleaned breath hose screen.

20 May 16
Date

D. Shanaver
Field Inspector's Signature

23 MAY 16 ✓
Date

Charles E. [Signature]
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:13

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 14214080A6
Standard Cyl #? 51
Standard Expiration? 07/05/2016
Oper No? 777777

Flow Cal. Date: 01/28/2010
Slope 738
Intercept -777574

IR Calibration Date: 01/28/2010

	3um	9um
0th Coef(*100):	-23128	-33609
1st Coef(*100):	272485	131724
2nd Coef(*100):	1490	1056
H2O adj(mg/l*10k):	432	614

***** Printer Test End *****


Operator Signature
DEB SHANAVER

Remarks:

Print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:32

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:33
02 Std. Sol.	0.017	13:35
03 Room Air	0.000	13:35
04 Std. Sol.	0.017	13:36
05 Room Air	0.000	13:36
06 Std. Sol.	0.018	13:37
07 Room Air	0.000	13:38

08 Sim Temp = 34.0°C

Simul Ser No = DR3374
Std Sol No = 201501B
County = 08 Oper No. = 777777



Operator Signature
DEB SHANAVER

Remarks:

low alcohol 0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:39
02 Std. Sol.	0.295	13:40
03 Room Air	0.000	13:40
04 Std. Sol.	0.295	13:41
05 Room Air	0.000	13:42
06 Std. Sol.	0.295	13:42
07 Room Air	0.000	13:43

08 Sim Temp = 34.0°C

Simul Ser No = DR3841
Std Sol No = 201603A
County = 08 Oper No. = 777777

Deb Shanaver
Operator Signature
DEB SHANAVER

Remarks:

High alcohol 0.300 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:45

Test	AC	Time
01 Room Air	0.000	13:46
02 *Subject Test	INT*	13:46
03 Room Air	0.000	13:47

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 14214080A6
Cyl No = 51
Expiration Date = 07/05/2016
County = 08 Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Interferent detected 0.10AC + 0.05% acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:48

Test	AC	Time
01 Diagnostic	OK	13:48
02 Room Air	RFI*	13:49
03 Room Air	0.000	13:49

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 14214080A6
Cyl No = 51
Expiration Date = 07/05/2016
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB SHANAVER

Remarks:

RFI test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:17
02 Std. Gas	0.079	13:17
03 Room Air	0.000	13:18
04 Std. Gas	0.079	13:18
05 Room Air	0.000	13:19
06 Std. Gas	0.080	13:19
07 Room Air	0.000	13:20

Lot No = 14214080A6
Cyl No = 51
Exp Date = 07/05/2016
County = 08 Oper No. = 777777

D.b Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:21
02 Std. Gas	0.080	13:21
03 Room Air	0.000	13:22
04 Std. Gas	0.079	13:22
05 Room Air	0.000	13:23
06 Std. Gas	0.080	13:23
07 Room Air	0.000	13:24

Lot No = 14214080A6
Cyl No = 51
Exp Date = 07/05/2016
County = 08 Oper No. = 777777



Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
05/20/2016 13:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:29
02 Std. Gas	0.079	13:29
03 Room Air	0.000	13:30
04 Std. Gas	0.080	13:30
05 Room Air	0.000	13:31
06 Std. Gas	0.079	13:31
07 Room Air	0.000	13:32

Lot No = 14214080A6
Cyl No = 51
Exp Date = 07/05/2016
County = 08 Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000