

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004203

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u> </u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5114</u> Lot No. <u> </u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u> </u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>709352</u>	Cyl No. <u>44</u>	Exp. Date <u>12/23/2011</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.079</u>	Test 5 <u>0.079</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.078</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 234 psi Reg 245 psi

M. Remarks/Maintenance record: annual inspection

28 May 2010
Date

Deb Kashner
Field Inspector's Signature

14 Jun 2010
Date

Mary Pearson
Reviewed by

80-004203
P. 57 of 65
KTH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 15:43

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 709352
Standard Cyl #? 44
Standard Expiration? 12/23/2011
Oper No? 777777

Flow Cal. Date: 01/28/2010
Slope 719
Intercept -711903

IR Calibration Date: 01/28/2010

	3um	9um
0th Coef(*100):	-17469	-28710
1st Coef(*100):	264216	135022
2nd Coef(*100):	2778	1081
H2O adj(mg/l*10k):	430	605

***** Printer Test End *****

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Print Test

Form 106-I8000

80-004203
P. 58 of 65
JAL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/27/2010 14:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:11
02 Std. Sol.	0.018	14:12
03 Room Air	0.000	14:13
04 Std. Sol.	0.018	14:13
05 Room Air	0.000	14:14
06 Std. Sol.	0.018	14:15
07 Room Air	0.000	14:15

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021 AC

Form 106-I8000

80-004203
P. 59 of 65
JK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/27/2010 14:04

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:04
02 Std. Sol.	0.293	14:05
03 Room Air	0.000	14:06
04 Std. Sol.	0.293	14:06
05 Room Air	0.000	14:07
06 Std. Sol.	0.293	14:08
07 Room Air	0.000	14:08

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Uncertainty test 0.298 AC

Form 106-18000

80-004203
P. 60 of 65
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 15:48

Test	AC	Time
01 Room Air	0.000	15:49
02 *Subject Test	INT*	15:49
03 Room Air	0.000	15:50

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 44
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

*Interferent check 0.100 AC
 + acetone*

Form 106-18000

80-004203
P. 61 of 65
JKM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 15:50

Test	AC	Time
01 Diagnostic	OK	15:51
02 Room Air	RFI*	15:51
03 Room Air	0.000	15:52

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 44
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-004203
P. 62 of 65
JML

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 13:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:56
02 Std. Gas	0.079	13:56
03 Room Air	0.000	13:56
04 Std. Gas	0.079	13:57
05 Room Air	0.000	13:57
06 Std. Gas	0.079	13:58
07 Room Air	0.000	13:58

Lot No = 709352

Cyl No = 44

Exp Date = 12/23/2011

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AE

Form 106-I8000

80-004203
P. 43465
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 14:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:21
02 Std. Gas	0.079	14:21
03 Room Air	0.000	14:22
04 Std. Gas	0.079	14:22
05 Room Air	0.000	14:22
06 Std. Gas	0.078	14:23
07 Room Air	0.000	14:23

Lot No = 709352
Cyl No = 44
Exp Date = 12/23/2011
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-18000

80-004203
p. 64 of 65
JK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = TOXL 8164.12.00 06/09
05/28/2010 14:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:42
02 Std. Gas	0.080	14:42
03 Room Air	0.000	14:43
04 Std. Gas	0.079	14:43
05 Room Air	0.000	14:44
06 Std. Gas	0.079	14:44
07 Room Air	0.000	14:45

Lot No = 709352

Cyl No = 44

Exp Date = 12/23/2011

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004203
p. 65 of 65
JK