



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Ray Weber</u>		
Location: <u>TRAIL CO. S.O.</u>		
Intoxilyzer® Serial Number: <u>80-004200</u>		
Gas Lot Number: <u>659358</u>	Gas Cylinder Number: <u>#45</u>	Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: <u>0.079</u> %		
Test 2: <u>0.079</u> %		
Test 3: <u>0.078</u> %		
(ATTACH TEST RECORD)		

Ray Weber
Chemical Test Operator Signature

5/29/2010
Date

Deb Kashur
Reviewed By

03 June 2010
Date

80-004200
P. 85 of 86
JRM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/29/2010 00:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	00:55
02 Std. Gas	0.079	00:55
03 Room Air	0.000	00:56
04 Std. Gas	0.079	00:56
05 Room Air	0.000	00:57
06 Std. Gas	0.078	00:57
07 Room Air	0.000	00:57

Lot No = 659358
Cyl No = 45
Exp Date = 09/03/2011
County = 49 Oper No. = 058405


Operator Signature
RAY G. WEBER

Remarks:

Form 106-I8000

ACA TEST

80-004200
P. 86 of 86
FM