



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Sgt Brian Crissler</i>		
Location <i>Belcourt Police Dept</i>	Intoxilyzer® Serial Number <i>80-003067</i>	
Gas Lot Number <i>068150</i>	Gas Cylinder Number <i>77</i>	Gas Expiration Date <i>9-11-2011</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. *081* AC
 - B. 0. *081* AC
 - C. 0. ~~080~~ AC *0.081 DK per BC old June 11*
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Sgt Brian Crissler</i>	Date <i>5-31-2011</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>09 June 2011</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = BELC 8164.13.00 06/09
05/31/2011 18:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:03
02 Std. Gas	0.081	18:03
03 Room Air	0.000	18:03
04 Std. Gas	0.081	18:04
05 Room Air	0.000	18:04
06 Std. Gas	0.081	18:05
07 Room Air	0.000	18:05

Lot No = 668180
Cyl No = 77
Exp Date = 09/21/2011
County = 63 Oper No. = 119802



Operator Signature
BRIAN F. CRISSLER

Remarks: *Change Cylinder out*

Form 106-I8000

ACA

80-003067
P. 100 of 100
VEM