



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Eugene M Malbert</i>		
Location <i>Adams Co SO HETT</i>	Intoxilyzer® Serial Number <i>80-004204</i>	
Gas Lot Number <i>05911080A1</i>	Gas Cylinder Number <i>7</i>	Gas Expiration Date <i>2-1-13</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.086* AC
 - B. *0.086* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Eugene M Malbert</i>	Date <i>5-31-11</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>18 July 2011</i>

*80-004204
P 112 of 113
JAM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004204
Location = HETT 8164.13.00 06/09
05/31/2011 19:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:42
02 Std. Gas	0.080	19:43
03 Room Air	0.000	19:43
04 Std. Gas	0.080	19:44
05 Room Air	0.000	19:44
06 Std. Gas	0.080	19:44
07 Room Air	0.000	19:45

Lot No = 05911080A1
Cyl No = 7
Exp Date = 02/01/2013
County = 01 Oper No. = 049130


Operator Signature
EUGENE M. MOLBERT

Remarks:

Form 106-I8000

80-004204
p. 113 of 113
JAN