



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Margy Pearson</u>		
Location: <u>Cass Co Jail</u>		
Intoxilyzer® Serial Number: <u>80-004205</u>		
Gas Lot Number: <u>709352</u>	Gas Cylinder Number: <u>14</u>	Expiration Date: <u>12/23/2011</u>
<input checked="" type="checkbox"/> New Cylinder Installed	<input type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%) Test 1: <u>0.081</u> % Test 2: <u>0.081</u> % Test 3: <u>0.080</u> % <p style="text-align: right;">(ATTACH TEST RECORD)</p>		

Margy Pearson
Chemical Test Operator Signature

3 Jan 2010
Date

Deb Kasher
Reviewed By

14 June 2010
Date

80-001205
P. 70 of 71
KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004205
Location = CASS 8164.13.00 06/09
06/03/2010 08:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:00
02 Std. Gas	0.081	08:01
03 Room Air	0.000	08:01
04 Std. Gas	0.081	08:02
05 Room Air	0.000	08:02
06 Std. Gas	0.080	08:02
07 Room Air	0.000	08:03

Lot No = 709352
Cyl No = 14
Exp Date = 12/23/2011
County = 09 Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Cyl Change ACA

Form 106-I8000

80-004205
P. 71 of 71
KAN