



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>CHARLES E. EDER (ANNUAL INSPECTION)</i>		
Location <i>GFSO (CORRECTIONAL)</i>	Intoxilyzer® Serial Number <i>80-003065</i>	
Gas Lot Number <i>33913080A4</i>	Gas Cylinder Number <i>047</i>	Gas Expiration Date <i>01/01/2016</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.080* AC
 - B. *0.080* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Charles E. Eder</i>	Date <i>03 JUN 14</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>04 June 14</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003065
Location = GFSO 8164.13.00 06/09
06/03/2014 10:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:26
02 Std. Gas	0.080	10:26
03 Room Air	0.000	10:27
04 Std. Gas	0.080	10:27
05 Room Air	0.000	10:28
06 Std. Gas	0.080	10:28
07 Room Air	0.000	10:29

Lot No = 33913080A4
Cyl No = 47
Exp Date = 01/01/2016
County = 18 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CYLINDER INSTALL*