



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Jamie M Huschka</u>		
Location: <u>Williams Co LFC</u>		
Intoxilyzer® Serial Number: <u>80-004189</u>		
Gas Lot Number: <u>659357</u>	Gas Cylinder Number: <u>45</u>	Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>76</u> %		
Test 2: 0. <u>76</u> %		
Test 3: 0. <u>76</u> %		
(ATTACH TEST RECORD)		

Jamie M Huschka
Chemical Test Operator Signature

16 June 2010
Date

Deb Kashner
Reviewed By

17 June 2010
Date

80-004189
p. 85 of 86
JMH


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = WILL 8164.13.00 06/09
06/16/2010 10:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:20
02 Std. Gas	0.076	10:21
03 Room Air	0.000	10:21
04 Std. Gas	0.076	10:22
05 Room Air	0.000	10:22
06 Std. Gas	0.076	10:22
07 Room Air	0.000	10:23

Lot No = 659357
Cyl No = 45
Exp Date = 09/03/2011
County = 53 Oper No. = 100406


Operator Signature
JAMIE M. HUSCHKA

Remarks:

Form 106-I8000

80-004189
P. 86 of 86
FM