



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>David J. Erickson</u>		
Location: <u>GRAF</u>		
Intoxilyzer® Serial Number: <u>80-003072</u>		
Gas Lot Number: <u>659358</u>	Gas Cylinder Number: <u>16</u>	Expiration Date: <u>09-03-11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>081</u> %		
Test 2: 0. <u>081</u> %		
Test 3: 0. <u>080</u> %		
(ATTACH TEST RECORD)		

David J. Erickson
Chemical Test Operator Signature

06-27-2010
Date

Deb Kashur
Reviewed By

29 June 2010
Date

80-003072
p. 89 of 90
JK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003072
Location = GRAF 8164.13.00 06/09
06/27/2010 23:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	23:46
02 Std. Gas	0.081	23:47
03 Room Air	0.000	23:47
04 Std. Gas	0.081	23:48
05 Room Air	0.000	23:48
06 Std. Gas	0.080	23:49
07 Room Air	0.000	23:49

Lot No = 659358
Cyl No = 16
Exp Date = 09/03/2011
County = 50 Oper No. = 119606


Operator Signature
DAVID J. ERICKSON

Remarks:

Form 106-I8000

80-003072
P. 90 of 90
JAN