



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Ray Weber</u>		
Location: <u>TRAIL CO S.O.</u>		
Intoxilyzer® Serial Number: <u>80-004200</u>		
Gas Lot Number: <u>659358</u>	Gas Cylinder Number: <u>#45</u>	Expiration Date: <u>09-03-11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: <u>0.080</u> %		
Test 2: <u>0.079</u> %		
Test 3: <u>0.078</u> %		
(ATTACH TEST RECORD)		

Ray Weber
 Chemical Test Operator Signature

7/1/10
 Date

Deb Kashur
 Reviewed By

16 July 2010
 Date

80-004200
 P. 87-88
 RM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
07/01/2010 21:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:46
02 Std. Gas	0.080	21:46
03 Room Air	0.000	21:47
04 Std. Gas	0.079	21:47
05 Room Air	0.000	21:48
06 Std. Gas	0.078	21:48
07 Room Air	0.000	21:49

Lot No = 659358
Cyl No = 45
Exp Date = 09/03/2011
County = 49 Oper No. = 058405


Operator Signature
RAY G. WEBER

Remarks:

Form 106-I8000

ACK TEST

80-004200
P. 88-88 pm