



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>MATHEW W. ANDERSON</i>		
Location <i>RICHLAND COUNTY JAIL</i>	Intoxilyzer® Serial Number <i>80-004206</i>	
Gas Lot Number <i>23411080A1</i>	Gas Cylinder Number <i>58</i>	Gas Expiration Date <i>10-01-13</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.080AC
 - B. 0.080AC
 - C. 0.080AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Sgt Mathew Anderson</i>	Date <i>7/05/12</i>
Reviewed By (Crime Laboratory Use Only) <i>Debi Kashner</i>	Date <i>13 July 2012</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
07/05/2012 17:18

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:18
02 Std. Gas	0.080	17:19
03 Room Air	0.000	17:19
04 Std. Gas	0.080	17:20
05 Room Air	0.000	17:20
06 Std. Gas	0.080	17:21
07 Room Air	0.000	17:21

Lot No = 23411080A1
Cyl No = 58
Exp Date = 10/01/2013
County = 39 Oper No. = 020401



Operator Signature
MATTHEW W. ANDERSON

Remarks: