



**ETHANOL BREATH STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Eugene M Molbert</i>		
Location <i>Adams Co SO HETT</i>	Intoxilyzer® Serial Number <i>80-003057</i>	
Gas Lot Number <i>116/1080 A2</i>	Gas Cylinder Number <i>112</i>	Gas Expiration Date <i>06-01-13</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0.080 AC  
 B. 0.080 AC  
 C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Eugene M Molbert</i>	Date <i>7-13-11</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>18 July 2011</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = HETT      8164.13.00 06/09  
07/13/2011      09:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:54
02 Std. Gas	0.080	09:54
03 Room Air	0.000	09:55
04 Std. Gas	0.080	09:55
05 Room Air	0.000	09:55
06 Std. Gas	0.080	09:56
07 Room Air	0.000	09:56

Lot No = 11611080A2  
Cyl No = 112  
Exp Date = 06/01/2013  
County = 01      Oper No. = 049130

  
Operator Signature  
EUGENE M. MOLBERT

Remarks:

Form 106-I8000

80-003057  
P. 149 of 149  
EM