



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) Pat Renz		
Location BSPD	Intoxilyzer® Serial Number <i>80-005947</i>	
Gas Lot Number 09216080A3	Gas Cylinder Number 35	Gas Expiration Date 05/05/2018

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 019 AC
 - B. 0. 078 AC
 - C. 0. 078 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Pat Renz</i>	Date <i>7-14-16</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date <i>25 July 16</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005947
Location = BSPD 8164.13.00 06/09
07/14/2016 16:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:34
02 Std. Gas	0.079	16:34
03 Room Air	0.000	16:34
04 Std. Gas	0.078	16:35
05 Room Air	0.000	16:35
06 Std. Gas	0.078	16:36
07 Room Air	0.000	16:36

Lot No = 09216080A3

Cyl No = 35

Exp Date = 05/05/2018

County = 08

Oper No. = 029710

Pat Renz

Operator Signature
PATRICK RENZ

Remarks: *New Cylinder Install*
New Test Record

Form 106-I8000