



ETHANOL GAS STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) JACQUIE DEMARCE		
Location FTOR	Intoxilyzer® Serial Number 80-005941	
Gas Lot Number 09216080#3	Gas Cylinder Number 26	Gas Expiration Date 5/5/2018

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.082 AC
 - B. 0.081 AC
 - C. 0.082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Jacqui Demarce</i>	Date 7-17-16 ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date 26 July 16


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005941
Location = FTOT 8164.13.00 06/09
07/17/2016 10:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:40
02 Std. Gas	0.082	10:40
03 Room Air	0.000	10:41
04 Std. Gas	0.081	10:41
05 Room Air	0.000	10:42
06 Std. Gas	0.082	10:42
07 Room Air	0.000	10:42

Lot No = 09216080A3
Cyl No = 26
Exp Date = 05/05/2018
County = 03 Oper No. = 100103


Operator Signature
JACQUIE DEMARCE

Remarks: