



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) CURTIS E. OLSON		
Location MCLEAN SO/GARRISON	Intoxilyzer® Serial Number 80-003072	
Gas Lot Number 80-003072 33913080 A4	Gas Cylinder Number 24	Gas Expiration Date 01-01-2016

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 081 AC
 - B. 0. 080 AC
 - C. 0. 081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Curtis E. Olson</i>	Date 7-19-2014
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date 21 July 14

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003072
Location = GARR 8164.13.00 06/09
07/19/2014 05:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	05:51
02 Std. Gas	0.081	05:51
03 Room Air	0.000	05:51
04 Std. Gas	0.080	05:52
05 Room Air	0.000	05:52
06 Std. Gas	0.081	05:53
07 Room Air	0.000	05:53

Lot No = 33913080A4
Cyl No = 24
Exp Date = 01/01/2016
County = 28 Oper No. = 059008



Operator Signature
CURTIS OLSON

Remarks: