

BREATH ALCOHOL INSTRUMENT CHAIN OF CUSTODY WORKSHEET

Instrument: 8000 Serial Number: 80-004184

Transported Out of the Crime Laboratory:

Agency Being Assigned Instrument: _____

Name/Agency of Person Transporting Instrument:

_____ (Print Name) _____ (Print Agency)

Signature of Person Receiving Instrument: _____

Date Receiving Instrument: ____/____/____ From: _____

Or Shipped Via:

USPS _____ UPS _____ FedEx Other _____

Date Sending Instrument: 21 July 2010 Air Bill No. 8697-8802-6296

Complete Air Bill (Identify Instrument/Serial Number) and attach a copy to this worksheet.

Transported to the Crime Laboratory:

Agency Returning Instrument: _____

Name/Agency of Person Returning Instrument:

_____ (Print Name) _____ (Print Agency)

Signature of Person Returning Instrument: _____

Date Returning Instrument: ____/____/____ Received By: _____

Or Received Via:

USPS _____ UPS _____ FedEx _____ Other _____

Date Receiving Instrument: ____/____/____ Air Bill No. _____

Attach shipping paperwork (if any) to this worksheet.

Instrument Transfer Information Entered Into the Equipment Database:

Date of Action 21 July 2010 By: Deb Kashner

80-004184
P. 75481
JH



Service Evaluation Form

This form **MUST** be completed and enclosed with instrument to be serviced.
Failure to complete and return this form may cause delays in service.

(Note: please ship items in their original shipping container(s) or a similar protective box.)

1. Contact information: Customer Number 58503 (contact Customer Service)
 Name Deb Kashur Phone: (701) 328-6159
 Fax: (701) 328-6185 Email: dkashur@nd.gov

2. Bill to Address: under warranty Ship to Address: ND Crime Laboratory
2641 E. Main Ave.
Bismarck, ND 58501-5044

3. Serial Number: 80-004184 Instrument Model: Intox 8000

4. Detailed Description of Problem:
In the field and the lab this intox doesn't always clear of alcohol w/higher AC's (0.290
"Purge Fails" are common. I've included a few tests. Also, it seems to only give
Purge Fails in CMS mode.
*****Hazardous Material Warning! - DO NOT return gas cylinders with instrument!*****

I Authorize Repairs Up To: All \$250 \$500 \$750 Other \$ under warranty
 Purchase Order Number (attach a copy of P.O. if applicable) _____

Authorized By:
Deb Kashur Forensic Scientist
 Name (Please Print) Title
Deb Kashur 21 July 2010
 Signature Date

Ship to:
 CMI, Inc.
 Attn: Service Dept.
 316 East Ninth Street
 Owensboro, KY 42303

No, please send estimate before repairs are made.
 Note: An estimate will be faxed before performing any repairs and may cause delays in service.
 An evaluation fee (\$79.00 or actual costs) will apply to estimates that are not repaired.

80-00 4184 5:09
 P. 7/21/10
 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = TOXL 8164.13.00 06/09
~~07/16/2010~~ 10:49

Test	AC	Time
01 Diagnostic	OK	10:50
02 Room Air	0.000	10:51
03 Subject Test 1	0.292	10:51
04 Room Air	PUR*	10:53
05 Room Air	0.000	10:53

lab test

*Purge Fail

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 47
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature
DEB KASHUR

Remarks:

Form 106-18000

*80-004184
P. 77 of 81 pp.
①*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = TOXL 8164.13.00 06/09
07/16/2010 11:03

Test	AC	Time
01 Diagnostic	OK	11:04
02 Room Air	0.000	11:04
03 Subject Test 1	0.289	11:05
04 Room Air	0.000 <i>too low</i>	11:07
05 Std. Gas	0.075	11:08
06 Room Air	0.000 <i>normally</i>	11:10
07 Subject Test 2	0.292 <i>0.079</i>	11:10
08 Room Air	0.000	11:11
09 Reported AC	0.289	11:05

Difference OK
No RFI Detected

Before this test I purged it (2A) for 1-2 min.

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 47
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature
DEB KASHUR

Remarks:

Form 106-I8000

*80-004184
P. 78 of 81
10/16
(2)*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = TOXL 8164.13.00 06/09
07/20/2010 07:56

Test	AC	Time
01 Diagnostic	OK	07:57
02 Room Air	0.000	07:58
03 Subject Test 1	0.292	07:58
04 Room Air	PUR*	08:00
05 Room Air	0.000	08:00

lab test

*Purge Fail

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 47
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature
DEB KASHUR

Remarks:

Form 106-I8000

80-004184
P. 79 of 81
132 (3)

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = TOXL 8164.13.00 06/09
07/20/2010 12:19

Test	AC	Time
01 Diagnostic	OK	12:20
02 Room Air	0.000	12:20
03 Subject Test 1	0.290	12:20
04 Room Air	0.000	12:23
05 Std. Gas	0.079 —	12:24
06 Room Air	0.000	12:25
07 Subject Test 2	0.290	12:26
08 Room Air	0.000	12:27
09 Reported AC	0.290	12:26

good lab test

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 47
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature
DEB KASHUR

Remarks:

Form 106-18000

*80-004184
P. 80481
KIL*

(4)

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = TOXL 8164.13.00 06/09
07/20/2010 12:32

Test	AC	Time
01 Diagnostic	OK	12:33
02 Room Air	0.000	12:33
03 Subject Test 1	0.288	12:34
04 Room Air	0.000	12:36
05 Std. Gas	0.077	12:37
06 Room Air	0.000	12:39
07 Subject Test 2	0.287	12:39
08 Room Air	0.000	12:40
09 Reported AC	0.287	12:39

good

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 709352
Cyl No = 47
Expiration Date = 12/23/2011
County = 08 Oper No. = 777777

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature
DEB KASHUR

Remarks:

Form 106-I8000

80-004184
P. 8/14/11
FILE
(5)