

BREATH ALCOHOL INSTRUMENT CHAIN OF CUSTODY WORKSHEET

Instrument: 8000 Serial Number: 80-004948

Transported Out of the Crime Laboratory:

Agency Being Assigned Instrument: _____

Name/Agency of Person Transporting Instrument:

(Print Name) _____ (Print Agency)

Signature of Person Receiving Instrument: _____

Date Receiving Instrument: ____/____/____ From: _____
Or Shipped Via:
USPS ____ UPS ____ FedEx ____ Other _____

Date Sending Instrument: ____/____/____ Air Bill No. _____

Complete Air Bill (Identify Instrument/Serial Number) and attach a copy to this worksheet.

Transported to the Crime Laboratory:

Agency Returning Instrument: _____

Name/Agency of Person Returning Instrument:

(Print Name) _____ (Print Agency)

Signature of Person Returning Instrument: _____

Date Returning Instrument: ____/____/____ Received By: _____
Or Received Via:
USPS ____ UPS ____ FedEx Other _____

Date Receiving Instrument: 23 July 2010 Air Bill No. no paperwork

Attach shipping paperwork (if any) to this worksheet.

Instrument Transfer Information Entered Into the Equipment Database:

Date of Action 23 July 2010 By: Deb Kashner

80-004948
P. 10/12
JRM

INTOX
PBT
TOXITEST

DATE
7-1-10
RECEIVED VIA
FELIX

CMI INC.

316 E. 9th St., Owensboro, KY 42303
270-685-6200 — Fax 270-685-6288
1-800-835-0690

SERVICE WORK ORDER

324330

B
I
L
L
T
O

ATTN:

CITY STATE ZIP

CUSTOMER NO. 585 052

PHONE

S
H
I
P
T
O

Office of Attorney General

ATTN: Crime Lab Division

2644 E. Main Ave.

Margy Pearson

CITY STATE ZIP
Bismarck ND 58501

PHONE
701-328-6159

MODEL #	S/N	DATE RCVD.	P.O. #	REQ.	BILL CODE	SERVICE LOCATION CODE
8000	80-004948				3	

PROD. CODE	EXTRA PARTS RCVD.
	BT, gas box, battery

DESCRIPTION OF PROBLEM
 "TEMP Prob - won't come out of 'Hot Ready'"

QUOTE
 WARRANTY
 SERVICE DEPT. NEW

- PRELIMINARY TEST DISMANTLE, INSPECT, DIAGNOSTIC TEST
- POWER SUPPLY TEST, INSPECT REPLACE
- MAIN PWB's - REMOVE & CLEAN, REINSTALL, TEST
- PWB's DIAGNOSE, REPAIR REPLACE
- PWB's - SURFACE CLEAN IN PLACE
- FILTER WHEEL - REMOVE, CLEAN, REINSTALL, CHECK CELL INTERIOR
- CHOPPER WHEEL MOTOR TEST & INSPECT REPLACE
- CELL - REMOVE, CLEAN, REINSTALL, REPLACE LENSES/GASKETS REQ'D
- IR LIGHT SOURCE REMOVE, REALIGN REPLACE
- PRINTER ASSY. DISASSEMBLE CLEAN & LUBRICATE
- REMOVE & REPLACE BELT REALIGNMENT
- REPLACE TEST PRINTER FUNCTIONS
- BATTERIES - REPLACE CPU RAM
- REPLACE PARTS REASSEMBLE CALIBRATION FINAL TEST
- MISC. _____
- PBT PRELIMINARY TESTING - PRESSURE, BATTERY, ZERO SET
- DISASSEMBLE REPLACE PARTS REASSEMBLE
- TEST CALIBRATE FINAL TEST
- TOXITEST INCOMING LEAK TEST VISUAL/MECH. INSPECTION
- TEST THERMOMETER & THERMOSTAT DISASSEMBLE
- REPLACE PARTS REASSEMBLE OUTGOING & LEAK TEST
- FINAL TEST
- MISC. _____

PARTS USED			
QTY	ID NO.	DESCRIPTION	CODE
1	021386FL	End block Assy	
1	021307	Breath hose Assy	

COMMENTS:
 I replaced the defective breath hose. In addition, I replaced the endblock Assy. I had this unit Re-card & final tested.

POA

HOURS: 4 hrs

TECH. DATE 7/13/10 CAL TECH. DATE 7/15/10

ID: #03471

SHIPPING DATE/CHARGES
 80-004948
 P 11/12

Certificate of Calibration

Intoxilyzer® 8000

*This is to certify that the calibration of **Intoxilyzer**® serial number 80-004948, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with the National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.*

Date 7/15/10

Signed Kristi Robinson
Technician

CMI INC.

80-004948
P. 12 of 12
vsm

316 East 9th Street
Owensboro, KY 42303
USA

Part No. 650519 CMI 7/9/09