



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) CHRISTOPHER A. SCHAEFER		
Location ROLLA	Intoxilyzer® Serial Number 86-004954	
Gas Lot Number 33913080A4	Gas Cylinder Number 39	Gas Expiration Date 01/01/2016

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0. 080 AC
- B. 0. 081 AC
- C. 0. 080 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature C.A. Schaefer	Date 07/25/2014
Reviewed By (Crime Laboratory Use Only) Deb Kashner	Date 28 July 2014

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = ROLL 8164.13.00 06/09
07/25/2014 17:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:50
02 Std. Gas	0.080	17:51
03 Room Air	0.000	17:51
04 Std. Gas	0.081	17:52
05 Room Air	0.000	17:52
06 Std. Gas	0.080	17:52
07 Room Air	0.000	17:53


Lot No = 33913080A4

Cyl No = 39

Exp Date = 01/01/2016

County = 40

Oper No. = 129973



Operator Signature
CHRISTOPHER SCHAEFER

Remarks: **NEW CYLINDER INSTALL**