



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) STOTLER, John W		
Location GFAF	Intoxilyzer® Serial Number 80-005954	
Gas Lot Number 13314080A3	Gas Cylinder Number 28	Gas Expiration Date 06/05/2016

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0.081 AC  
 B. 0.081 AC  
 C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>John W Stotler</i>	Date 07/28/2014
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date 28 July 14

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005954  
Location = GFAP      8164.13.00 06/09  
07/28/2014      14:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:21
02 Std. Gas	0.081	14:21
03 Room Air	0.000	14:22
04 Std. Gas	0.081	14:22
05 Room Air	0.000	14:23
06 Std. Gas	0.081	14:23
07 Room Air	0.000	14:24

Lot No = 13314080A3

Cyl No = 28

Exp Date = 06/05/2016

County = 65

Oper No. = 031121



Operator Signature  
JOHN STOTLER

Remarks:

Form 106-18000