



# ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>FAVE L. AXT</i>		
Location <i>HARVEY PD</i>	Intoxilyzer® Serial Number <i>80-002666</i>	
Gas Lot Number <i>33913080 A4</i>	Gas Cylinder Number <i>43</i>	Gas Expiration Date <i>01/01/2016</i>

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0.079 AC

B. 0.079 AC

C. 0.079 AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <i>Fave Axt</i>	Date <i>7-31-14</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>04 Aug 2014</i>

THIS REPORT IS THE PROPERTY OF THE STATE OF NORTH CAROLINA. IT IS TO BE RETURNED TO THE OFFICE OF ATTORNEY GENERAL, CRIME LABORATORY DIVISION, TOXICOLOGY SECTION, BREATH ALCOHOL PROGRAM, 100 SOUTH BRIDGE STREET, RALEIGH, NC 27601. (919) 733-2200


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002666  
Location = HARV      8164.13.00 06/09  
07/31/2014      23:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	23:28
02 Std. Gas	0.079	23:28
03 Room Air	0.000	23:29
04 Std. Gas	0.079	23:29
05 Room Air	0.000	23:30
06 Std. Gas	0.079	23:30
07 Room Air	0.000	23:31

Lot No = 33913080A4  
Cyl No = 43  
Exp Date = 01/01/2016  
County = 52      Oper No. = 110201

  
\_\_\_\_\_  
Operator Signature  
FAYE AXT

Remarks: