



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) Collins, Travis L.		
Location HETT	Intoxilyzer® Serial Number 80-004952	
Gas Lot Number 33913080A4	Gas Cylinder Number 13	Gas Expiration Date 01/01/2016

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 078 AC
 - B. 0. 078 AC
 - C. 0. 078 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>T Collins</i>	Date 08/07/2014
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date 13 Aug 2014

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004952
Location = HETT 8164.13.00 06/09
08/07/2014 07:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	07:34
02 Std. Gas	0.078	07:34
03 Room Air	0.000	07:35
04 Std. Gas	0.078	07:35
05 Room Air	0.000	07:36
06 Std. Gas	0.078	07:36
07 Room Air	0.000	07:37

Lot No = 33913080A4
Cyl No = 13
Exp Date = 01/01/2016
County = 01 Oper No. = 020533



Operator Signature
TRAVIS COLLINS

Remarks:

Form 106-I8000