OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Program
Breath Alcohol Instrument Chain of Custody Worksheet

BREATH ALCOHOL INSTRUMENT CHAIN OF CUSTODY WORKSHEET

Instrument: <u> </u>	tox 8000		Serial Number: _	80-004	184
Transported Out of the	Crime Laboratory:		VALUE	· · · · · · · · · · · · · · · · · · ·	
Agency Being Assigned Instrument:					
Name/Agency of Perso	n Transporting Insti	rument:			
(Print Name)					_ (Print Agency)
Signature of Person Receiving Instrument:					
Date Receiving Instrum	ent:/	/	From:	***************************************	
Or Shipped Via:					
USPS	UPS Fe	edEx	Other		
Date Sending Instrume	nt:/		Air Bill No		
Complete Air Bill (Identify Instrument/Serial Number) and attach a copy to this worksheet.					
Transported to the Crim	e Laboratory:				
Agency Returning Instrument:					
Name/Agency of Person Returning Instrument:					
	(Prin	nt Name) _		(1	Print Agency)
Signature of Person Returning Instrument:					
Date Returning Instrument:/Received By:					
Or Received Via:					
USPS UPS FedEx _ / Other					
Date Receiving Instrument: 09 1 Aug 2010 Air Bill No. no paperwork					
Attach shipping paperwork (if any) to this worksheet.					
Instrument Transfer Information Entered Into the Equipment Database:					
Date of Action <u>091 Aug 1 2010</u> By: <u>Deb Kashur</u>					

86-004/84 P. 83-485

Worksheet: BrW-012

SERVICE WORK ORDER INTOX 316 E. 9th St., Owensboro, KY 42303 PBT 270-685-6200 — Fax 270-685-6288 RECEIVED VIA TOXITEST 1-800-835-0690 В S ATTN: CITY T CUSTOMER NO. PHONE MODEL# DATE ROVD. P.O. # REQ. BILL CODE SERVICE LOCATION CODE PROD. CODE DESCRIPTION OF PROBLEM WARRANTY **□** QUOTE ☐ SERVICE DEPT. ☐ NEW Pulge Fails PRELIMINARY TEST DISMANTLE, INSPECT, DIAGNOSTIC TEST PARTS USED POWER SUPPLY ☐ TEST, INSPECT ☐ REPLACE DESCRIPTION QTY ID NO. CODE MAIN PWB's - REMOVE & CLEAN, REINSTALL, TEST 020982 Sunzer MSSY REPLACE ☐ DIAGNOSE, REPAIR □ PWB's - SURFACE CLEAN IN PLACE ☐ FILTER WHEEL - REMOVE, CLEAN, REINSTALL, CHECK CELL INTERIOR CHOPPER WHEEL MOTOR ☐ TEST & INSPECT ☐ REPLACE CELL - REMOVE, CLEAN, REINSTALL, REPLACE LENSES/GASKETS REQ'D IR LIGHT SOURCE ☐ REMOVE, REALIGN —☐ REPLACE PRINTER ASSY. DISASSEMBLE CLEAN & LUBRICATE ☐ REMOVE & REPLACE BELT ☐ REALIGNMENT ☐ REPLACE ☐ TEST PRINTER FUNCTIONS BATTERIES - REPLACE ☐ RAM. □ CPU REPLACE PARTS PREASSEMBLE ☐ FINAL TEST ☐ MISC. PBT | PRELIMINARY TESTING - PRESSURE, BATTERY, ZERO SET ☐ DISASSEMBLE ☐ REPLACE PARTS REASSEMBLE ☐ TEST ☐ CALIBRATE ☐ FINAL TEST TOXITEST I INCOMING LEAK TEST VISUAL/MECH. INSPECTION ☐ TEST THERMOMETER & THERMOSTAT ☐ DISASSEMBLE ☐ REPLACE PARTS ☐ REASSEMBLE ☐ OUTGOING & LEAK TEST ☐ FINAL TEST ☐ MISC. COMMENTS: Arthurd Exhaust Block Rullices Parts listed SHIPPING DATE/CHARGES 5-21 # 53552 80-004184 HOURS: 4.0HA

8/2/15

FORM#

CMI 8WO 1 REV (-)

SPECIAL SHIPPING INSTRUCTIONS

COD

☐ AMOUNT \$

Certificate of Calibration Intoxityzer® 8000

This is to certify that the calibration of Intoxilyzer® serial number RO-OO4184, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with the National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

Date 8/3/10 Signed Knisti Robinson



80-004184 f.85.4.85

316 East 9th Street Owensboro, KY 42303 USA

Part No. 650519 CMI 7/9/09