



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (5-2011)

Chemical Test Operator Name (Print): <u>Misialek, Kyle D.</u>		
Location: <u>(Cros) Divide Co. S.O.</u>	Intoxilyzer® Serial Number: <u>80-004184</u>	
Gas Lot Number: <u>13314080AD</u>	Gas Cylinder Number: <u>57</u>	Gas Expiration Date: <u>06/05/16</u>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information. (Level 1, Function S)
- 2. Perform an ACA Test. (Level 1, Function C)

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081):

- A. 0.079 AC
- B. 0.078 AC
- C. 0.079 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN 50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN 50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN 59282, Form 105-G).
 - B. ACA Test Record.

[Signature]
 Chemical Test Operator Signature

8/12/14
 Date

Deb Kashner
 Reviewed By (Crime Laboratory Use Only)

14 Aug 14
 Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = CROS 8164.13.00 06/09
08/12/2014 10:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:43
02 Std. Gas	0.079	10:43
03 Room Air	0.000	10:44
04 Std. Gas	0.078	10:44
05 Room Air	0.000	10:44
06 Std. Gas	0.079	10:45
07 Room Air	0.000	10:45

Lot No = 13314080A2
Cyl No = 57
Exp Date = 06/05/2016
County = 12 Oper No. = 130670



Operator Signature
KYLE MISIALEK

Remarks: Test OK