



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Tyler Wolfe</i>		
Location <i>BOTT</i>	Intoxilyzer® Serial Number <i>80-004189</i>	
Gas Lot Number <i>13314080A2</i>	Gas Cylinder Number <i>52</i>	Gas Expiration Date <i>6/05/2016</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0.079 AC  
 B. 0.080 AC  
 C. 0.079 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Tyler Wolfe</i>	Date <i>8-25-14</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Debs Kashu</i>	Date <i>05 Sept 14</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004189  
Location = BOTT                 8164.13.00 06/09  
08/25/2014                        13:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:18
02 Std. Gas	0.079	13:18
03 Room Air	0.000	13:18
04 Std. Gas	0.080	13:19
05 Room Air	0.000	13:19
06 Std. Gas	0.079	13:20
07 Room Air	0.000	13:20

Lot No = 13314080A2  
Cyl No = 52  
Exp Date = 06/05/2016  
County = 05                            Oper No. = 060901

  
\_\_\_\_\_  
Operator Signature  
TYLER WOLFE

Remarks: