



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Darcy W Aberle</i>		
Location <i>Williston LEC</i>	Intoxilyzer® Serial Number <i>80-004950</i>	
Gas Lot Number <i>29614080A2</i>	Gas Cylinder Number <i>9</i>	Gas Expiration Date <i>12/05/16</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. 084 AC  
 B. 0. 084 AC  
 C. 0. 084 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Darcy Aberle</i>	Date <i>09-06-15</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>09 Sept 15</i>

*Aberle found the AC result to be reading high, yet still in range, so he kept running tests to see what the results would be. There are 4 tests in this group.  
 DK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004950  
Location = WILL      8164.13.00 06/09  
09/06/2015      14:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:25
02 Std. Gas	0.084	14:25
03 Room Air	0.000	14:26
04 Std. Gas	0.084	14:26
05 Room Air	0.000	14:26
06 Std. Gas	0.084	14:27
07 Room Air	0.000	14:27

Lot No = 29614080A2  
Cyl No = 9  
Exp Date = 12/05/2016  
County = 53      Oper No. = 010101

  
Operator Signature  
DARCY ABERLE

Remarks:



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Darcy Aberle</i>		
Location <i>Williston LEC</i>	Intoxilyzer® Serial Number <i>80-004950</i>	
Gas Lot Number <i>29614080A2</i>	Gas Cylinder Number <i>9</i>	Gas Expiration Date <i>12/05/16</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. 084 AC  
 B. 0. 084 AC  
 C. 0. 084 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Darcy Aberle</i>	Date <i>09/06/15</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>09 Sept 15</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004950  
Location = WILL      8164.13.00 06/09  
09/06/2015      14:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:30
02 Std. Gas	0.084	14:30
03 Room Air	0.000	14:31
04 Std. Gas	0.084	14:31
05 Room Air	0.000	14:31
06 Std. Gas	0.084	14:32
07 Room Air	0.000	14:32

Lot No = 29614080A2  
Cyl No = 9  
Exp Date = 12/05/2016  
County = 53      Oper No. = 010101

  
Operator Signature  
DARCY ABERLE

Remarks:



**ETHANOL GAS STANDARD CYLINDER REPORT**  
OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Darcy Aberle</i>		
Location <i>Williston LEC</i>	Intoxilyzer® Serial Number <i>80-004950</i>	
Gas Lot Number <i>13314080A2</i>	Gas Cylinder Number <i>51</i>	Gas Expiration Date <i>06/05/16</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.084*AC
  - B. *0.084*AC
  - C. *0.083*AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Darcy Aberle</i>	Date <i>09/06/15</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date <i>09 Sept 15</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004950  
Location = WILL      8164.13.00 06/09  
09/06/2015      14:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:35
02 Std. Gas	0.084	14:35
03 Room Air	0.000	14:36
04 Std. Gas	0.084	14:36
05 Room Air	0.000	14:37
06 Std. Gas	0.083	14:37
07 Room Air	0.000	14:38

Lot No = 13314080A2

Cyl No = 51

Exp Date = 06/05/2016

County = 53

Oper No. = 010101



Operator Signature  
DARCY ABERLE

Remarks:

Form 106-I8000



**ETHANOL GAS STANDARD CYLINDER REPORT**  
OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Darcy Aberle</i>		
Location <i>Williston LEC</i>	Intoxilyzer® Serial Number <i>80-004950</i>	
Gas Lot Number <i>29614080A2</i>	Gas Cylinder Number <i>9</i>	Gas Expiration Date <i>12/05/16</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.085* AC
  - B. *0.084* AC
  - C. *0.085* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Darcy Aberle</i>	Date <i>09/06/15</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>09 Sept 15</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004950  
Location = WILL                8164.13.00 06/09  
09/06/2015                      14:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:40
02 Std. Gas	0.085	14:40
03 Room Air	0.000	14:40
04 Std. Gas	0.084	14:41
05 Room Air	0.000	14:41
06 Std. Gas	0.085	14:42
07 Room Air	0.000	14:42

Lot No = 29614080A2  
Cyl No = 9  
Exp Date = 12/05/2016  
County = 53                      Oper No. = 010101



Operator Signature  
DARCY ABERLE

Remarks: