



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>David Felond</u>		
Location: <u>Kidder SO.</u>		
Intoxilyzer® Serial Number: <u>80-004186</u>		
Gas Lot Number: <u>659357</u>	Gas Cylinder Number: <u>37</u>	Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: <u>0.081</u> %		
Test 2: <u>0.080</u> %		
Test 3: <u>0.080</u> %		
(ATTACH TEST RECORD)		

David Felond
Chemical Test Operator Signature

15 Sep 10
Date

Deb Kashner
Reviewed By

16 Sept 10
Date

80-004186
P. 89 of 90
JKM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = STEE 8164.13.00 06/09
09/15/2010 14:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:05
02 Std. Gas	0.081	14:05
03 Room Air	0.000	14:06
04 Std. Gas	0.080	14:06
05 Room Air	0.000	14:07
06 Std. Gas	0.080	14:07
07 Room Air	0.000	14:08

Lot No = 659357
Cyl No = 37
Exp Date = 09/03/2011
County = 22 Oper No. = 119404



Operator Signature
DAVID R. FELAND

Remarks:

Form 106-I8000

80-004186
P. 90 of 90
FM