



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Lt. Tracy L. Meidinger #258</u>		
Location: <u>UND-PD</u>		
Intoxilyzer® Serial Number: <u>80-004934</u>		
Gas Lot Number: <u>659358</u>	Gas Cylinder Number: <u>31</u>	Expiration Date: <u>09-30-11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>080</u> %		
Test 2: 0. <u>079</u> %		
Test 3: 0. <u>080</u> %		
(ATTACH TEST RECORD)		

Lt. Tracy L. Meidinger #258
 Chemical Test Operator Signature

10-7-10
 Date

Deb Kashur
 Reviewed By

18 Oct 2010
 Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = UNDP 8164.13.00 06/09
10/07/2010 22:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:11
02 Std. Gas	0.080	22:12
03 Room Air	0.000	22:12
04 Std. Gas	0.079	22:13
05 Room Air	0.000	22:13
06 Std. Gas	0.080	22:13
07 Room Air	0.000	22:14

Lot No = 659358
Cyl No = 31
Exp Date = 09/03/2011
County = 18 Oper No. = 119409

Tracy L. Meidinger #258

Operator Signature
TRACY L. MEIDINGER

Remarks:

Form 106-I8000

ACA October 2010

80-004934
p. 58 of 58
FM