



# ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) CHRISTOPHER SCHAEFER		
Location ROLLA	Intoxilyzer® Serial Number 80-004954	
Gas Lot Number 16512080 A Z	Gas Cylinder Number 050	Gas Expiration Date 08/01/2014

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.083 AC
  - B. 0.082 AC
  - C. 0.083 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date 10/24/2013 ✓
Reviewed By (Crime Laboratory Use Only) Deb Kashun	Date 01 NOV 2013

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004954  
Location = ROLL      8164.13.00 06/09  
10/24/2013      10:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:38
02 Std. Gas	0.083	10:38
03 Room Air	0.000	10:38
04 Std. Gas	0.082	10:39
05 Room Air	0.000	10:39
06 Std. Gas	0.083	10:40
07 Room Air	0.000	10:40

Lot No = 16512080A2

Cyl No = 050

Exp Date = 08/01/2014

County = 40

Oper No. = 129973



Operator Signature  
CHRISTOPHER SCHAEFER

Remarks:

Form 106-I8000