



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Ray Weller</u>		
Location: <u>TRAIL CO. SO.</u>		
Intoxilyzer® Serial Number: <u>80-004200</u>		
Gas Lot Number: <u>1541008041</u>	Gas Cylinder Number: <u>#15</u>	Expiration Date: <u>09/01/2012</u>
<input checked="" type="checkbox"/> New Cylinder Installed <input type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%) Test 1: 0. <u>076</u> % Test 2: 0. <u>076</u> % Test 3: 0. <u>076</u> %		
(ATTACH TEST RECORD)		

Ray Weller
 Chemical Test Operator Signature

10/25/2010
 Date

Deb Kashner
 Reviewed By

27 Oct 2010
 Date

80-004200
 P. 93 of 94
 jcm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
10/25/2010 15:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:24
02 Std. Gas	0.076	15:24
03 Room Air	0.000	15:25
04 Std. Gas	0.076	15:25
05 Room Air	0.000	15:26
06 Std. Gas	0.076	15:26
07 Room Air	0.000	15:27

Lot No = 15410080A1

Cyl No = 15

Exp Date = 07/01/2012

County = 49

Oper No. = 058405



Operator Signature
RAY G. WEBER

Remarks:

Form 106-I8000

80-004200
p. 94 of 94
JMN