



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Jeremy Gravis</i>		
Location <i>Burke County Sheriff's Office</i>		Intoxilyzer® Serial Number <i>80-004940</i>
Gas Lot Number <i>16512080A2</i>	Gas Cylinder Number <i>57</i>	Gas Expiration Date <i>8-1-2014</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.078 AC
 - B. 0.078 AC
 - C. 0.077 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Jeremy Gravis</i>	Date <i>10/26/2013</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Keshur</i>	Date <i>01 NOV 2013</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004940
Location = BURK 8164.13.00 06/09
10/26/2013 17:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:41
02 Std. Gas	0.078	17:41
03 Room Air	0.000	17:42
04 Std. Gas	0.078	17:42
05 Room Air	0.000	17:42
06 Std. Gas	0.077	17:43
07 Room Air	0.000	17:43

Lot No = 16512080A2

Cyl No = 57

Exp Date = 08/01/2014

County = 07

Oper No. = 130230


Operator Signature
JEREMY GROHS

Remarks:

Cylinder Change

Form 106-I8000