



OFFICE OF ATTORNEY GENERAL
 Crime Laboratory Division
 2641 East Main Avenue
 P.O. Box 937
 Bismarck, North Dakota 58502-0937

TEL 701-328-6159

ETHANOL BREATH STANDARD CYLINDER INSTALLATION FORM

(Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.)

Chemical Test Operator Name: Ray Weller

Location: TRAIL CO. S.O.

Intoxilyzer Serial No.: 80-004191

Lot Number: 659358

Cylinder Number: #45

Expiration Date: 09/03/2011

Instrument Test Results (Report to 3 Digits; Example, 0.081%)

Test 1: 0. 080 %

Test 2: 0. 081 %

Test 3: 0. 082 %

(ATTACH TEST RECORD)

Ray Weller
 Chemical Test Operator Signature
Deb Kashur

11-07-2009
 Date Cylinder Changed
17 Nov 2009

Form 105-G
 Toxicology Laboratory
 Crime Laboratory Division

80-004191
 P-1 of 2 ym

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = HILL 8164.12.00 06/09
11/07/2009 03:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	03:24
02 Std. Gas	0.080	03:24
03 Room Air	0.000	03:25
04 Std. Gas	0.081	03:25
05 Room Air	0.000	03:26
06 Std. Gas	0.082	03:26
07 Room Air	0.000	03:26

Lot No = 659358

Cyl No = 45

Exp Date = 09/03/2011

County = 49

Oper No. = 058405



Operator Signature
RAY G. WEBER

Remarks:

ACA TEST

Form 106-18000

80-004191
P. 2 of 2 jwr