

INTOXILYZER® 5000 G AND 8000
INSTALLATION AND REPAIR CHECKOUT

Serial No.: 80 103066 Instrument Location: Dickinson ZSC ✓

Reason for Install/Repair: _____ Install After Lab Repair/Inspection
x Other (specify) New Machine From Toxicology

Check When Done:

- A. Surge Protector Installed/Properly Grounded.
- B. Telephone Line Connected to Intoxilyzer® KB.
- C. Breath Tube Heated.
- D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- E. Ethanol Breath Standard: Connected _____ Installed.
(Include Form 105-G if New Cylinder Installed.)
- F. Review/Enter Gas Setup (Level 1, Function S). (Intoxilyzer® 8000 Only.)
- G. Print Test (Level 1, Function P). (Attach Test Record.)
- H. ACA Test (Level 1, Function C). (Attach Test Record.)
- I. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- J. Repair and/or Maintenance Performed: _____

11-18-09
Date
24 NOV 2009
Date

[Signature]
Field Inspector's Signature
Deb Kashur
Reviewed By

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = DKPD 8164.13.00 06/09
11/18/2009 21:28

Test	AC	Time
01 Diagnostic	OK	21:30
02 Room Air	RFI*	21:30
03 Room Air	0.000	21:31

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/08/1987
Sub Sex = Male Weight = 150
Test = OTH Cit = 1
Dr. Lic. = ND/DIS871456
Lot No = 659358
Cyl No = 23
Expiration Date = 09/03/2011
County = 45 Oper No. = 068332

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DAVID M. WALLACE

Remarks:

Install

Form 106-I8000

80-003066
P. 2 of 4
FW

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = DKPD 8164.13.00 06/09
11/18/2009 21:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:21
02 Std. Gas	0.081	21:21
03 Room Air	0.000	21:21
04 Std. Gas	0.081	21:22
05 Room Air	0.000	21:22
06 Std. Gas	0.082	21:23
07 Room Air	0.000	21:23

Lot No = 659358

Cyl No = 23

Exp Date = 09/03/2011

County = 45

Oper No. = 068332



Operator Signature
DAVID M. WALLACE

Remarks:

Install

Form 106-I8000

80-003066
P. 3 of 4
JW

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = DKPD 8164.13.00 06/09
11/18/2009 21:19

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 23
Standard Expiration? 09/03/2011
Oper No? 068332

Flow Cal. Date: 04/08/2008
Slope 696
Intercept -643215

IR Calibration Date: 06/04/2008
 3um 9um

0th Coef(*100): -11861 -9861
1st Coef(*100): 275638 138672
2nd Coef(*100): 165 672
H2O adj(mg/l*10k): 357 250

***** Printer Test End *****



Operator Signature
DAVID M. WALLACE

Remarks: *Install*