OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Program
Breath Alcohol Instrument Chain of Custody Worksheet

BREATH ALCOHOL INSTRUMENT CHAIN OF CUSTODY WORKSHEET

Instrument: 8000	Serial Number: _	80-003057
Transported Out of the Crime Laboratory:		
Agency Being Assigned Instrument:		
Name/Agency of Person Transporting Instrument:		
(Print Name	e)	(Print Agency)
Signature of Person Receiving Instrument:		
Date Receiving Instrument://	From:	
Or Shipped Via:		
USPS UPS FedEx _	Other	
Date Sending Instrument://	Air Bill No	
Complete Air Bill (Identify Instrument/Serial Number) and attach a copy to this worksheet.		
Transported to the Crime Laboratory:		
Agency Returning Instrument:		
Name/Agency of Person Returning Instrument:		
(Print Name)		(Print Agency)
Signature of Person Returning Instrument:		
Date Returning Instrument://	Received By:	
Or Red	ceived Via:	
USPS UPS FedEx _\	Other	
Date Receiving Instrument: 13 1 Dec 1 2011 Air Bill No.		
Attach shipping paperwork (if any) to this worksheet.		
Instrument Transfer Information Entered Into the Equipment Database:		
Date of Action 13 / Dec 1 20(1 By: Deb Kashur		

Worksheet: BrW-012

Certificate of Calibration *Intoxityzer* 8000

This is to certify that the calibration of INIONITYZET® serial number 80-03057, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with the National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

Date 11-29-10

Signed

Technician



80-003057 P.118 4/19 316 East 9th Street Owensboro, KY 42303 USA

Part No. 650519 CMi 7/9/0:

SERVICE WORK ORDER NTOX ulesiul 316 E. 9th St., Owensboro, KY 42303 PBT 335354 270-685-6200 — Fax 270-685-6288 RECEIVED VIA TOXITEST 1-800-835-0690 F/-LUV& Attorney General В ATTN: S ATTN: Home CITY STATE Т T CUSTOMER NO. PHONE MODEL# DATE ROVD. P.O. # REQ. BILL CODE SERVICE LOCATION CODE 8000 WARRANTY QUOTE SERVICE DEPT. I NEW ☐ DISMANTLE, INSPECT, DIAGNOSTIC TEST ☑ PRELIMINARY TEST **PARTS USED** POWER SUPPLY ☐ TEST, INSPECT REPLACE QTY ID NO. DESCRIPTION CODE ☐ MAIN PWB's - REMOVE & CLEAN, REINSTALL, TEST W1388 Recirculation persons PWB's DIAGNOSE, REPAIR REPLACE ☐ PWB's - SURFACE CLEAN IN PLACE ☐ FILTER WHEEL - REMOVE, CLEAN, REINSTALL, CHECK CELL INTERIOR CHOPPER WHEEL MOTOR TEST & INSPECT REPLACE ☐ CELL - REMOVE, CLEAN, REINSTALL, REPLACE LENSES/GASKETS REQ'D IR LIGHT SOURCE ☐ REMOVE, REALIGN REPLACE PRINTER ASSY. ☐ DISASSEMBLE ☐ CLEAN & LUBRICATE ☐ REALIGNMENT ☐ REMOVE & REPLACE BELT ☐ REPLACE ☐ TEST PRINTER FUNCTIONS BATTERIES - REPLACE ☐ CPU ☐ RAM FINAL TEST ☐ REASSEMBLE ☐ CALIBRATION ☐ REPLACE PARTS ☐ MISC. PBT PRELIMINARY TESTING - PRESSURE, BATTERY, ZERO SET ☐ DISASSEMBLE ☐ REPLACE PARTS ☐ REASSEMBLE ☐ TEST ☐ CALIBRATE ☐ FINAL TEST TOXITEST ☐ INCOMING LEAK TEST ☐ VISUAL/MECH. INSPECTION ☐ TEST THERMOMETER & THERMOSTAT ☐ DISASSEMBLE ☐ REPLACE PARTS ☐ REASSEMBLE ☐ OUTGOING & LEAK TEST ☐ FINAL TEST ☐ MISC. COMMENTS: I Uprished the Problem & I settled the both a simulator colonnet & Re-set the all goins. SHIPPING DATE/CHARGES

CMI 8WO 1 REV (-) FORM#

80-003057 P.119 4/19

Mh

SPECIAL SHIPPING INSTRUCTIONS

☐ AMOUNT \$